



Product Datasheet

TNFRSF25 antibody (orb542066)


Description	Rabbit polyclonal antibody to TNFRSF25
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human TNFRSF25(NP_683866.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	60kDa
Uniprot ID	Q93038
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of rat breast tissue.



Immunofluorescence analysis of human colon tissue.



Western blot analysis of extracts of various tissues.