



Product Datasheet

TRIM15 antibody (orb542130)


Description	Rabbit polyclonal antibody to TRIM15
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TRIM15(NP_150232.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	52kDa
Uniprot ID	Q9C019
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human uterine tissue.



Immunohistochemical staining of human venous tissue.



Immunohistochemical staining of rat liver tissue.