



Product Datasheet

TRIM27 antibody (orb542134)


Description	Rabbit polyclonal antibody to TRIM27
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TRIM27 (NP_006501.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	59kDa
Uniprot ID	P14373
Entrez	5987
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:10-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...



Immunohistochemical staining of mouse st...



Immunohistochemical staining of human es...