


Product Datasheet

TRIM45 antibody (orb542140)

Description	Rabbit polyclonal antibody to TRIM45
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human TRIM45(NP_001139107.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	64KDa
Uniprot ID	Q9H8W5
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...