


Product Datasheet

TSG101 antibody (orb542165)

Description	Rabbit polyclonal antibody to TSG101
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	A synthetic peptide of human TSG101/VPS23(NP_006283.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	31kDa/43kDa
Uniprot ID	Q99816
Dilution Range	WB: 1:200-1:2000, IP: 1:20-1:50
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...