



Product Datasheet

Tubulin beta 3 antibody (orb542183)


Description	Rabbit Polyclonal antibody to Tubulin beta 3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human β III-Tubulin/ β 3-Tubulin/ β III-Tubulin/ β 3-Tubulin.
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	55kDa
Uniprot ID	Q13509
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200, IP: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-251MG C...



Immunofluorescence analysis of rat brain...



Immunofluorescence analysis of mouse bra...