



Product Datasheet

TXNRD1 antibody (orb542199)


Description	Rabbit polyclonal antibody to TXNRD1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human Thioredoxin reductase 1 (TXNRD1)(Q16881).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65kDa
Uniprot ID	Q16881
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse testis



Immunohistochemical staining of human placenta



Western blot analysis of extracts of various tissues