

Product Datasheet

Human LSP1 protein (orb54220)

Description	Recombinant human LSP1 protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal GST-tagged
Note	For research use only
Application notes	N-terminal GST-tagged: N-terminal GST-tagged1-198AA: 1-339AFull Length : Full Length of BC001785
Protein Sequence	MAEASSDPGAEEEREELLGPTAQWSVEDEEEAVHEQCQHERDRQLQAQDEEGG GHVPERPKQEMLLSLKPSEAPELDEDEGFGDWSQRPEQRQQHEGAQGTLD SG EPPQCRSPEGEQEDRPGLHAYEKEDSDEVHLEELSLSKEGPGPEDTVQDNLGAA GAEEEEQEEHQKCQQRTPSPLVLEGTIEQSSPPLSPTTKLIDRTESLNRSIEKSNS VKKSQPDLPISKIDQWLEQYTQAIETAGRTPKLARQASIELPSMAVASTKSRWET GEVQAQSAAKTPSCKDIVAGDMSKSLWEQKGGSKTSSTIKSTPSGKRYKFVAT GHGKYEKVLVEGGPAP
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	64.2 kDa
Uniprot ID	P33241
Expiration Date	6 months from date of receipt.