

Product Datasheet

UBE2M antibody (orb542215)

Description	Rabbit polyclonal antibody to UBE2M
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human UBE2M(NP_003960.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	21kDa
Uniprot ID	P61081
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.