



## Product Datasheet

### UQCR10 antibody (orb542257)


|                            |                                                          |
|----------------------------|----------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to UQCR10                     |
| <b>Reactivity</b>          | Human, Mouse, Rat                                        |
| <b>Conjugation</b>         | Unconjugated                                             |
| <b>Tested Applications</b> | ICC, IF, IHC, WB                                         |
| <b>Immunogen</b>           | Recombinant fusion protein of human UQCR10(NP_037519.2). |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.       |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.              |
| <b>Note</b>                | For research use only                                    |
| <b>Isotype</b>             | IgG                                                      |
| <b>Clonality</b>           | Polyclonal                                               |
| <b>MW</b>                  | 10kDa                                                    |
| <b>Uniprot ID</b>          | <b>Q9UDW1</b>                                            |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, IHC-P: 1:50-1:200                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                          |



Immunohistochemical staining of rat hear...



Immunohistochemical staining of rat ovar...



Western blot analysis of extracts of var...