



Product Datasheet

VCP antibody (orb542298)


Description	Rabbit polyclonal antibody to VCP
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human VCP(NP_009057.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	90KDa
Uniprot ID	P55072
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-2 OS cells...



Immunofluorescence analysis of C6 cells ...



Immunohistochemical staining of human esophagus...