



Product Datasheet

VEGF antibody (orb542305)


Description	Rabbit polyclonal antibody to VEGF
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human VEGFA(NP_001165099.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	23KDa,27KDa,45KDa
Uniprot ID	P15692
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HUVEC cells...



Western blot analysis of extracts of var...



Western blot analysis of extracts of var...