


Product Datasheet

YY1 antibody (orb542383)

Description	Rabbit polyclonal antibody to YY1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, WB
Immunogen	Recombinant fusion protein of human YY1(NP_003394.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65KDa
Uniprot ID	P25490
Dilution Range	WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...