


Product Datasheet

ZC3H12A antibody (orb542395)

Description	Rabbit polyclonal antibody to ZC3H12A
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human ZC3H12A(NP_079355.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	66kDa
Uniprot ID	Q5D1E8
Dilution Range	WB: 1:1000-1:3000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...