



Product Datasheet

ZNF581 antibody (orb542447)


Description	Rabbit polyclonal antibody to ZNF581
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human ZNF581(NP_057619.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22-24kDa
Uniprot ID	Q9P0T4
Dilution Range	WB: 1:1000-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-2 OS cells



Immunofluorescence analysis of NIH-3T3 cells



Immunofluorescence analysis of C6 cells