

Product Datasheet

ZSCAN4 antibody (orb542467)

Description	Rabbit polyclonal antibody to ZSCAN4
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human ZSCAN4(NP_689890.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	49kDa
Uniprot ID	Q8NAM6
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.