



Product Datasheet

STAT1 (Phospho-S727) antibody (orb542497)


Description	Rabbit polyclonal antibody to STAT1 (Phospho-S727)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic phosphorylated peptide around S727 of human STAT1 (NP_009330.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	87kDa
Uniprot ID	P42224
Entrez	6772
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100, IP: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat ovar...



Immunohistochemical staining of rat test...



Immunohistochemical staining of rat live...