



Product Datasheet

RIPK1 antibody (orb542509)


Description	Rabbit polyclonal antibody to RIPK1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	A synthetic peptide of human RIPK1/RIP(NP_003795.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	75KDa/38KDa
Uniprot ID	Q13546
Dilution Range	WB: 1:500-1:2000, IP: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse sp...



Immunohistochemical staining of rat lung...



Immunohistochemical staining of human to...