



Product Datasheet

IL1B antibody (orb542527)


Description	Rabbit polyclonal antibody to IL1B
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant protein of human IL1 β .
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	37kDa
Uniprot ID	P01584
Entrez	3553
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of rat spleen...



Immunofluorescence analysis of mouse spleen...



Western blot analysis of extracts of various cell lines...