

Product Datasheet

FGF1 antibody (orb542543)

Description	Rabbit polyclonal antibody to FGF1
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human FGF1(NP_000791.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	6kDa/17kDa
Uniprot ID	P05230
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.