

Product Datasheet

COL11A1 antibody (orb542565)

Description	Rabbit polyclonal antibody to COL11A1
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human COL11A1(NP_542196.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	167kDa/176kDa/181kDa/182kDa
Uniprot ID	P12107
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.