

Product Datasheet

AHI1 antibody (orb542631)

Description	Rabbit polyclonal antibody to AHI1
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human AHI1(NP_001128304.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70kDa/120kDa/137kDa
Uniprot ID	Q8N157
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.