

Product Datasheet

RNASE12 antibody (orb542668)

Description	Rabbit polyclonal antibody to RNASE12
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human RNASE12 (NP_001019993.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q5GAN4
Entrez	493901
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.