

Product Datasheet

SCGB2A1 antibody (orb542715)

Description	Rabbit polyclonal antibody to SCGB2A1
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SCGB2A1(NP_002398.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	10kDa
Uniprot ID	075556
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.