

## Product Datasheet

### Cadherin 5 antibody (orb542716)

|                            |                                                               |
|----------------------------|---------------------------------------------------------------|
| <b>Description</b>         | Rabbit Polyclonal antibody to Cadherin 5                      |
| <b>Reactivity</b>          | Human, Mouse                                                  |
| <b>Conjugation</b>         | Unconjugated                                                  |
| <b>Tested Applications</b> | IHC, WB                                                       |
| <b>Immunogen</b>           | Recombinant fusion protein of human VE Cadherin(NP_001786.2). |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal,50% glycerol,pH7.3.                 |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                   |
| <b>Note</b>                | For research use only                                         |
| <b>Isotype</b>             | IgG                                                           |
| <b>Clonality</b>           | Polyclonal                                                    |
| <b>MW</b>                  | 140KDa                                                        |
| <b>Uniprot ID</b>          | <b>P33151</b>                                                 |
| <b>Dilution Range</b>      | IF: 1:50-1:200                                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                               |