

Product Datasheet

SCG3 antibody (orb542734)

Description	Rabbit polyclonal antibody to SCG3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SCG3 (NP_037375.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q8WXD2
Entrez	29106
Dilution Range	WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.