

Product Datasheet

C1QTNF4 antibody (orb542736)

Description	Rabbit polyclonal antibody to C1QTNF4
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human C1QTNF4 (NP_114115.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q9BXJ3
Entrez	114900
Dilution Range	WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.