

Product Datasheet

LYPD1 antibody (orb542748)

Description	Rabbit polyclonal antibody to LYPD1
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human LYPD1(NP_653187.3).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	9kDa/15kDa/16kDa/17kDa
Uniprot ID	Q8N2G4
Dilution Range	WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.