

Product Datasheet

Lysozyme antibody (orb542752)

Description	Rabbit Polyclonal antibody to Lysozyme
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Lysozyme (Lysozyme (LYZ)) (NP_000230.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P61626
Entrez	4069
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.