

Product Datasheet

USP6 antibody (orb542765)

Description	Rabbit polyclonal antibody to USP6
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human USP6(NP_004496.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	89kDa/121kDa/158kDa
Uniprot ID	P35125
Dilution Range	WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.