

Product Datasheet

C1QA antibody (orb542777)

Description	Rabbit polyclonal antibody to C1QA
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human C1QA(NP_057075.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	26kDa
Uniprot ID	P02745
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.