

Product Datasheet

CPB2 antibody (orb542779)

Description	Rabbit polyclonal antibody to CPB2
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human CPB2(NP_001863.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40kDa/48kDa
Uniprot ID	Q96IY4
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.