

Product Datasheet

ASCC2 antibody (orb542790)

Description	Rabbit polyclonal antibody to ASCC2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human ASCC2 (NP_115580.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q9H1I8
Entrez	84164
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.