

Product Datasheet

SLC5A6 antibody (orb542800)

Description	Rabbit polyclonal antibody to SLC5A6
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Recombinant fusion protein of human SLC5A6 (NP_066918.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q9Y289
Entrez	8884
Dilution Range	IF: 1:50-1:100
Expiration Date	12 months from date of receipt.