

Product Datasheet

CBP80 antibody (orb542805)

Description	Rabbit polyclonal antibody to CBP80
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide of human CBP80(NP_002477.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	100KDa
Uniprot ID	Q09161
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.