

Product Datasheet

CGRP antibody (orb542806)

Description	Rabbit Polyclonal antibody to CGRP
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human CALCA(NP_001029124.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	13kDa/15kDa
Uniprot ID	P06881 , P01258
Dilution Range	IF: 1:50-1:200
Expiration Date	12 months from date of receipt.