

Product Datasheet

BLOC1S1 antibody (orb542812)

Description	Rabbit polyclonal antibody to BLOC1S1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human BLOC1S1 (NP_001478.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P78537
Entrez	2647
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.