

Product Datasheet

KIR3DL1 antibody (orb542817)

Description	Rabbit polyclonal antibody to KIR3DL1
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human KIR3DL1 (NP_037421.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	60kDa
Uniprot ID	P43629
Entrez	3811
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.