

Product Datasheet

RANTES antibody (orb542854)

Description	Rabbit Polyclonal antibody to RANTES
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human CCL5/RANTES (NP_002976.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	16kDa
Uniprot ID	P13501
Entrez	6352
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.