

Product Datasheet

TAF1B antibody (orb542857)

Description	Rabbit polyclonal antibody to TAF1B
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TAF1B(NP_005671.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	39kDa/60kDa/68kDa
Uniprot ID	Q53T94
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.