

Product Datasheet

ZNF169 antibody (orb542861)

Description	Rabbit polyclonal antibody to ZNF169
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human ZNF169 (NP_919301.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q14929
Entrez	169841
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.