

Product Datasheet

HSD17B6 antibody (orb542867)

Description	Rabbit polyclonal antibody to HSD17B6
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human HSD17B6(NP_003716.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36kDa
Uniprot ID	O14756
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.