

Product Datasheet

DAZAP2 antibody (orb542868)

Description	Rabbit polyclonal antibody to DAZAP2
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human DAZAP2(NP_055579.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	17kDa
Uniprot ID	Q15038
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.