

Product Datasheet

NECAP2 antibody (orb542875)

Description	Rabbit polyclonal antibody to NECAP2
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human NECAP2(NP_060560.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	28kDa
Uniprot ID	Q9NVZ3
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.