

Product Datasheet

RHBDD1 antibody (orb542877)

Description	Rabbit polyclonal antibody to RHBDD1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human RHBDD1(NP_115652.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	12kDa/35kDa
Uniprot ID	Q8TEB9
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.