

Product Datasheet

TAF6L antibody (orb542878)

Description	Rabbit polyclonal antibody to TAF6L
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant fusion protein of human TAF6L(NP_006464.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70kDa
Uniprot ID	Q9Y6J9
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.