

Product Datasheet

ZNF585A antibody (orb542892)

Description	Rabbit polyclonal antibody to ZNF585A
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human ZNF585A(NP_001275729.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	46kDa/87kDa
Uniprot ID	Q6P3V2
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.