

Product Datasheet

SUMF1 antibody (orb542896)

Description	Rabbit polyclonal antibody to SUMF1
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SUMF1(NP_877437.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	31kDa/37kDa/38kDa/40kDa/46kDa
Uniprot ID	Q8NBK3
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.