

## Product Datasheet

### DHRS11 antibody (orb542924)

|                            |                                                          |
|----------------------------|----------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to DHRS11                     |
| <b>Reactivity</b>          | Human, Mouse                                             |
| <b>Conjugation</b>         | Unconjugated                                             |
| <b>Tested Applications</b> | WB                                                       |
| <b>Immunogen</b>           | Recombinant fusion protein of human DHRS11(NP_077284.2). |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal,50% glycerol,pH7.3.            |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.              |
| <b>Note</b>                | For research use only                                    |
| <b>Isotype</b>             | IgG                                                      |
| <b>Clonality</b>           | Polyclonal                                               |
| <b>MW</b>                  | 20kDa/28kDa                                              |
| <b>Uniprot ID</b>          | <b>Q6UWP2</b>                                            |
| <b>Dilution Range</b>      | WB: 1:500-1:2000                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                          |