

Product Datasheet

ALDH3B2 antibody (orb542935)

Description	Rabbit polyclonal antibody to ALDH3B2
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human ALDH3B2(NP_000686.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48KDa
Uniprot ID	P48448
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.